


Hemo Control

How to perform a test



1 Ensure that the analyzer is ready for use. Refer to Hemo Control Operator's Manual for details.



2 Take one cuvette from the canister. Close the lid of the canister immediately. Only use the handle and do not touch the optical eye of the microcuvette.



3 Disinfect the puncture site and allow it to dry. Make the incision on the side of the finger tip. Support finger and apply light pressure while making the incision.



4 Apply only light pressure towards the fingertip and discard the first drops (1-3). Make sure there is free blood flow before filling the cuvette.



5 Fill the microcuvette in one step by holding the tip on to the middle of the blood drop. Ensure there are no air bubbles present.



6 Remove any blood from the outside of the microcuvette. Ensure that no blood is drawn out of the microcuvette during this process.



7 Open the cuvette holder by pulling it out carefully. Insert the cuvette with the raised wing pointing to the right of the cuvette holder.



8 The measurement will begin when the cuvette holder is closed.



9 The result is displayed on the screen after ca. 25 - 60 seconds. Press OK to leave the result screen.

Manufacturer

EKF diagnostic GmbH
Ebendorfer Chaussee 3
39179 Barleben, Germany
Tel +49 (0)39203 5110

Sales

EKF Diagnostics Holdings plc
Avon House, 19 Stanwell Road,
Penarth, Cardiff, UK, CF64 2EZ
Tel +44 (0)29 20 710570

EKF | Diagnostics
for life