

Lactate Scout 4



Sports
Performance



NEW

Connection
to heart rate
monitors

NEW

E-paper display for
easier viewing on
the move

NEW

Ergonomic design,
smaller, lighter,
more robust

Lactate Scout 4

Accurate lactate
measurement to drive
performance improvement

- ◆ Result in 10 seconds
- ◆ Only a drop of fingerprick blood
- ◆ Simple to use
- ◆ Smartphone and tablet compatible Bluetooth®
- ◆ Compensated for hematocrit

<p>Easy to use</p>	<ul style="list-style-type: none"> • Insert/remove sensor to turn on/off • Simplified navigation to enable flexibility whilst training • Only 0.2 µl of capillary blood required • Pre-calibrated sensors • Device calibrated by simple coding • Lactate measurement value, date/time, measuring mode, temperature, heart rate and memory ID on a single display • Automatic self test
<p>Fast and accurate</p>	<ul style="list-style-type: none"> • Enzymatic amperometric detection method • Results within 10 seconds • Measuring range: 0.5 - 25 mmol/L • Compensates for the influence of low and high hematocrit levels • Imprecision: <ul style="list-style-type: none"> - Hct-range 35 - 50%: 0.5 - 6.7 mmol/L blood lactate standard deviation ≤ 0.2 mmol/L, 6.8 - 25 mmol/L blood lactate CV $\leq 3\%$ - Hct range 20 - <35%, >50 - 70%: 0.5 - 7.5 mmol/L blood lactate standard deviation ≤ 0.3 mmol/L, 7.6 - 25 mmol/L blood lactate CV $\leq 4\%$ • Test solutions available for function control
<p>Performance measurement</p>	<ul style="list-style-type: none"> • Single and step-test measurements (resting/exercise/recreation) • Connection to heart rate monitors compatible to Bluetooth® Low Energy for linking to lactate value • Stopwatch function • Lactate Scout Assistant software available for performance management
<p>Practical and reliable</p>	<ul style="list-style-type: none"> • Stores up to 500 results • 1,000 tests using just 2 x CR2450 batteries • Integrated Bluetooth® Low Energy connectivity • Pocket-size: 91 mm (h) x 46 mm (w) x 21 mm (d) • Lightweight: 60 g • Operating range: 10 - 45°C and max. 85% humidity

Catalog No.	Name
7023-0441-0246	Lactate Scout 4 Solo
7023-0441-0568	Lactate Scout 4 Start Set



1 Remove sensor and place into analyser.



2 Prick finger and collect blood by touching with sensor.



3 Result appears in 10 seconds.

Distributed by

Manufacturer
SensLab GmbH
Bautzner Str. 67
04347 Leipzig
Germany

+49 (0) 39203 511 0
sales@ekfdiagnostics.com

ekfdiagnostics.com

EKF | Diagnostics
for life