


# SAFE AND SECURE SINGLE USE LANCETS FOR OBTAINING CAPILLARY BLOOD SAMPLES



Easy to use

Safe

Innovative

Efficient

**ibabyLance**



**EKF**  
**DIAGNOSTICS**

SCIENCE MADE SIMPLE

<b>EASY TO USE</b>	<ul style="list-style-type: none"> <li>• One step procedure</li> <li>• No arming required</li> </ul> <p><b>SurgiLance:</b></p> <ul style="list-style-type: none"> <li>• Activated by pressure on finger</li> </ul> <p><b>BabyLance:</b></p> <ul style="list-style-type: none"> <li>• Reduced trigger activation force, significantly reduced risk of bruising</li> </ul>
<b>SAFE</b>	<ul style="list-style-type: none"> <li>• Needle or blade safely retracted and concealed before and after use</li> <li>• Latex-free, hypoallergenic</li> <li>• Once used the lancets are rendered inoperative, providing added safety for patient and clinician</li> </ul>
<b>INNOVATIVE</b>	<ul style="list-style-type: none"> <li>• Compact, ergonomic design</li> <li>• Colour-coded for easy recognition of penetration depth and size</li> <li>• Optimal blood flow for different applications and skin types</li> </ul>
<b>EFFICIENT</b>	<ul style="list-style-type: none"> <li>• High-speed delivery and penetration method</li> <li>• Minimised patient discomfort</li> <li>• Comfortable to use</li> </ul>

Medipurpose Safety Lancets give Point-of-Care providers a quick and effective testing method, minimising pain while maximising patient protection.

The pre-armed needle or blade is concealed before and after use to prevent cross-contamination and needlestick injuries.

Their unmatched reliability makes Medipurpose Safety Lancets extremely cost-effective. Include lancets in your test kits. Or make them available with the user's preferred PoC equipment.

**SurgiLance and BabyLance: Safety made simple.**

## SurgiLance

Colour	Depth	Gauge	Blood flow
Grey	1.8mm	G23	10-20µl
Orange	2.2mm	G22	20-40µl
Pink	2.8mm	G21	40-60µl



## BabyLance

Colour	Depth	Code
Pink	0.85mm	BLP Preemie
Blue	1.0mm	BLP Newborn



Distributed by

**Manufacturer**  
**MediPurpose™ Private Limited**  
 3883 Rogers Bridge Rd. NW  
 Suite 501  
 Duluth, Georgia 30097  
**Tel.:** +1-770-448-9493

**Distributor**  
**EKF Diagnostics Holdings plc**  
 Avon House, 19 Stanwell Road,  
 Penarth, Cardiff, UK, CF64 2EZ  
**Tel.:** +44 (0)29 20 710570  
**www.ekfdiagnostics.com**

**EKF**  
**DIAGNOSTICS**  
 SCIENCE MADE SIMPLE