



## PrimeStore® MTM Sample Collection Kit

- Sample inactivation media designed for the safe transportation and handling of pathogenic samples
- Convenient CE marked sample collection kit designed for home use and mass on-site sampling
- Test for multiple nucleic acid species from one sample
- Sample stability for up to 7 days at 2-25 °C
- No cold chain required for sample collection or transportation

## Simple and safe sample collection

The PrimeStore® MTM sample collection kit enables convenient and safe collection, transportation and handling of pathogenic samples including Covid-19 and Influenza.

### Bulk sample collection kit

The PrimeStore® MTM bulk sample collection kits are available in a pack size of 75 kits. Each sample collection kit contains the following components:

- 1 x UN3373 specimen collection bag
- 1 x Cryovial containing 1.5mL of PrimeStore molecular transport media
- 1 x Sterile CE marked nasopharyngeal flocked swab
- 1 x Absorbent pad
- 1 x Instructions for use



Product code: PS-MTM-SK

### Postal sample collection kit

The PrimeStore® MTM postal sample collection kit is ideal for home sampling, the kit includes the following components:

- 1 x UN3373 specimen collection bag
- 1 x Cryovial containing 1.5mL of PrimeStore molecular transport media
- 1 x Sterile CE marked nasopharyngeal flocked swab
- 1 x Absorbent pad
- 1 x Instructions for use
- 1 x UN3373 postal box with security seal for sample return
- 1 x Postal pouch for collection kit delivery



Product code: PS-MTM-PSK

Kit components	Specification
PrimeStore MTM	1.5mL of PrimeStore MTM inactivation media in a 5.0mL self standing cryovial with external thread and leak-proof caps. 24 months shelf-life at 2-25 °C (36-77 °F) from date of manufacture
Flocked swab	Sterile, flocked standard tip nasopharyngeal swab with an 80mm breakpoint
Absorbent Pad	Absorbent Sheet 75x125mm capable of absorbing 50ml of liquid
Specimen collection bag	For transport of category B biological samples to UN3373 P650 packaging standards
UN3373 postal box with security seal for sample return	Sample transport box for mailing of category B biological samples to UN3373 P650 packaging standards
Postal pouch for collection kit delivery	Padded envelope to contain postal kit contents for delivery to subject's address

Distributed by

#### Sales

EKF Diagnostics

☎ 029 2071 0570

✉ sales@ekfdiagnostics.com

www.ekfdiagnostics.com

#### Manufacturer

Stanbio Laboratory

☎ +44 (0)29 20710570

✉ info@ekfdiagnostics.com

www.ekfdiagnostics.com



Diagnostics  
for life