

Quo-Test®



Quo-Test®

Easy and reliable HbA1c measurement

- ◆ Fully automated analyzer
- ◆ Fast and accurate
- ◆ Easy to use
- ◆ Results in just 4 minutes

Simple, accurate and reliable HbA1c measurement

<p>Easy to use</p>	<ul style="list-style-type: none"> • Intuitive, step by step on screen instructions minimise training time • Step by step instructions on display • Multi-lingual menu • Just 4µl blood required from finger prick or venous sample
<p>Efficient data handling</p>	<ul style="list-style-type: none"> • Barcode reader included. Only one scan per lot of cartridges required • Stores up to 7000 results for efficient patient data management • User selectable dual reporting. User can select % DCCT, IFCC mmol/mol, eAG mg/dl or eAG mmol/l. • USB Port • Optional label printer available
<p>Fast and accurate</p>	<ul style="list-style-type: none"> • Results within 4 minutes • Boronate affinity method which is recognized as interference free • Measuring range: 4–15% A1c DCCT • Imprecision: CV <3% at 7% A1c DCCT • Unaffected by Hb variants, which do not result in reduced erythrocyte life span • Traceable to the IFCC reference method
<p>Ergonomic design</p>	<ul style="list-style-type: none"> • Compact: 205 mm(h) x 205 mm(w) x 135 mm(d) • Lightweight: 1.3kg

Catalog No.	Name
3121-0011_VS	<ul style="list-style-type: none"> • Quo-Test Analyzer Pack



1 Collect 4 µl blood sample.



2 Insert sample collector into the Quo-Test cartridge.



3 Insert cartridge into analyser. Results displayed within 4 minutes.

Distributed by

Manufacturer

EKF-Diagnostic GmbH
 Ebendorfer Chaussee 3
 39179 Barleben
 Germany

+49 (0) 39203 511 0
 sales@ekfdiagnostics.com

www.ekfdiagnostics.com

EKF | Diagnostics
 for life