



QuStick™ Strep A Rapid Test

Early diagnosis of Strep A infection means more cost-effective treatment and increased patient satisfaction

- ◆ Enhanced sensitivity
- ◆ Superior culture classification
- ◆ Easy-to-use
- ◆ Results within five minutes



STANBIO®



Diagnostics
for life

QuStick™ Strep A Rapid Test

Rapid, qualitative detection of group A streptococcal antigen directly from patient throat specimen.

Enhanced sensitivity

97% sensitivity compared to culture

QuStick™ Strep A's ultra-sensitivity means that it requires fewer colonies to yield a positive result. Which means that even patients presenting less obvious symptoms can be detected

Superior culture classification (low counts)

91 % @ Rare (< 5 colonies)

100 % @ 1+ (< 10 colonies)

Easy-to-use

CLIA waived, requires minimal training

Fast

Results in five minutes

Economical. Convenient.

QuStick™ eliminates waste with room temperature storage and a two year shelf life from date of manufacture

Complete

The kit includes everything you need

Each kit includes:

- Individually pouched test strips
- Reagent A and reagent B
- Positive and negative external controls
- Sterile throat swabs
- Extraction tubes
- Instructions for use
- Handy pictorial quick guide showing step-by-step procedure



Reference Nos.
6000-025 (25 tests)
6000-050 (50 tests)

Distributed by



CODE*:
87880QW

ISO 13485:
2003 Certified

Stanbio
1261 North Main Street
Boerne, Texas
78006 USA

+1 830.249.0772
866.782.6246 (U.S. Toll Free)
sales2@ekfdiagnostics.com

www.ekfdiagnostics.com

STANBIO® **EKF** | Diagnostics for life

* Exemplar code provided. Refer to regional CMS fee schedule for relevant codes and reimbursement rates.