



## RELY® Mono

Detect Infectious Mononucleosis quickly and reliably with a proven, easy-to-use, rapid test

- ◆ Two step test provides results in five minutes
- ◆ Convenient storage temperature (2 - 30 °C)
- ◆ Uses whole blood, serum or plasma

# RELY® Mono

RELY® Mononucleosis Rapid Test provides a simple, fast and cost-effective test for Infectious Mononucleosis. Symptoms include fever, sore throat and swollen lymph glands and mimic other serious infections that need different and urgent treatment.

## Enhanced sensitivity\*

98%

## Specificity\*

>99%

## Reliable

Procedural and external controls provided

## Sample

Whole blood (fingerstick or venipuncture), serum or plasma

## Easy-to-use

Requires minimal training  
CLIA waived: whole blood

## Cost effective

Convenient room temperature/refrigerated storage and extended dating help manage inventory during seasonal shifts in testing

**Storage:** 2-30 °C (36-86 °F)

**Shelf life:** 18 months from date of manufacture

Quick and reliable detection of Mononucleosis means improved patient care

- Minimize use of unnecessary antibiotics
- Quick results allow for immediate therapy to ease discomfort and shorten patient's wait time
- Early IgM detection helps reduce the risk of infecting others



Reference No.  
6200-020 (20 tests)

Distributed by



CODE\*\*:  
86308QW (whole blood)  
86308 (serum/plasma)

ISO 13485: 2003 Certified

Stanbio  
1261 North Main Street  
Boerne, Texas  
78006 USA

+1 830.249.0772  
1-800-531-5535 (USA Toll Free)  
sales2@ekfdiagnostics.com

[www.ekfdiagnostics.com](http://www.ekfdiagnostics.com)

STANBIO®

EKF | Diagnostics  
for life

\*For all sample types  
\*\* Exemplar code provided. Refer to CMS fee schedule for relevant codes and reimbursement rates.