


SAFE AND SECURE SINGLE USE LANCETS FOR OBTAINING CAPILLARY BLOOD SAMPLES



Easy to use

Safe

Innovative

Efficient

ibabyLance



EKF
DIAGNOSTICS

SCIENCE MADE SIMPLE

EASY TO USE	<ul style="list-style-type: none"> • One step procedure • No arming required <p>SurgiLance:</p> <ul style="list-style-type: none"> • Activated by pressure on finger <p>BabyLance:</p> <ul style="list-style-type: none"> • Reduced trigger activation force, significantly reduced risk of bruising
SAFE	<ul style="list-style-type: none"> • Needle or blade safely retracted and concealed before and after use • Latex-free, hypoallergenic • Once used the lancets are rendered inoperative, providing added safety for patient and clinician
INNOVATIVE	<ul style="list-style-type: none"> • Compact, ergonomic design • Colour-coded for easy recognition of penetration depth and size • Optimal blood flow for different applications and skin types
EFFICIENT	<ul style="list-style-type: none"> • High-speed delivery and penetration method • Minimised patient discomfort • Comfortable to use

Medipurpose Safety Lancets give Point-of-Care providers a quick and effective testing method, minimising pain while maximising patient protection.

The pre-armed needle or blade is concealed before and after use to prevent cross-contamination and needlestick injuries.

Their unmatched reliability makes Medipurpose Safety Lancets extremely cost-effective. Include lancets in your test kits. Or make them available with the user's preferred PoC equipment.

SurgiLance and BabyLance: Safety made simple.

SurgiLance

Colour	Depth	Gauge	Blood flow
Grey	1.8mm	G23	10-20µl
Orange	2.2mm	G22	20-40µl
Pink	2.8mm	G21	40-60µl



BabyLance

Colour	Depth	Code
Pink	0.85mm	BLP Preemie
Blue	1.0mm	BLP Newborn



Distributed by

Manufacturer
MediPurpose™ Private Limited
 3883 Rogers Bridge Rd. NW
 Suite 501
 Duluth, Georgia 30097
Tel.: +1-770-448-9493

Distributor
EKF Diagnostics Holdings plc
 Avon House, 19 Stanwell Road,
 Penarth, Cardiff, UK, CF64 2EZ
Tel.: +44 (0)29 20 710570
www.ekfdiagnostics.com

EKF
DIAGNOSTICS
 SCIENCE MADE SIMPLE