

**Trade name:** Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution

**Current version :** 8.0.0, issued: 22.03.2013

**Replaced version:** 7.1.0, issued: 14.03.2013

**Region:** GB

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1 Product identifier

Trade name

**Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution**

**System solution**

**Order No. (Epp.-No.)**

**0201-0002-024**

**0201-0002-025**

**0201-0002-026**

**0201-0002-030 (0030 322.002)**

**Micro test tubes with hemolyzing solution**

**Order No. (Epp.-No.)**

**0209-0100-003**

**0209-0100-005 (6668 054.003)**

**0209-0100-011**

**0209-0100-012**

**0209-0100-013**

**0209-0100-014**

**0209-0100-050**

**0209-0100-122 (6660 050.009)**

**0209-0100-1244**

**0209-0100-135 (6660 040.003)**

**0209-0100-138 (6660 140.008)**

**0209-0100-1414**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses of the substance or mixture**

Glucose/Lactate hemolyzing solution:

Bottles and micro test tubes with hemolyzing solution with or without end-to-end capillaries are used to carry out 1 + 50 dilution [hemolysis] of blood for the determination of glucose and lactate. These 1 + 50 dilutions can be measured using the devices BIOSEN.

Glucose/Lactate system solution:

Glucose/Lactate system solution used in the devices BIOSEN to condition the biosensor and as a rinsing and transport solution. Glucose/Lactate system solution can also be used as a hemolyzing solution (sample dilution: 1 + 50) in the BIOSEN devices.

**Trade name:** Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution

**Current version :** 8.0.0, issued: 22.03.2013

**Replaced version:** 7.1.0, issued: 14.03.2013

**Region:** GB

## 1.3 Details of the supplier of the safety data sheet

### Address

EKF-diagnostic GmbH  
Ebendorfer Chaussee 3  
D-39179 Barleben

Telephone no. +49 (0)3 92 03 / 7 85 - 0  
Fax no. +49 (0)3 92 03 / 7 85 - 16

### Information provided by / telephone

EKF-diagnostic GmbH  
Telephone: +49 (0)39203 / 785 - 0 (7 am - 4 pm on workdays)

### Advice on Safety Data Sheet

support@ekf-diagnostic.de

## 1.4 Emergency telephone number

+49 (0) 551 / 19240 Poison Information Centre North(Giftinformationszentrale Nord)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

#### Classification information

This product does not meet the classification and labelling criteria given in the Dangerous Preparations Directive (1999/45/EC; DPD)

### 2.2 Label elements

Not relevant

### 2.3 Other hazards

No special hazards have to be mentioned. Please observe the information given in this safety data sheet.

## SECTION 3: Composition/information on ingredients

### 3.1 Substances

Not applicable. The product is not a substance.

### 3.2 Mixtures

#### Chemical characterization

aqueous solution

#### Hazardous ingredients

This product does not contain substances to be mentioned according to EU Regulation No 1907/2006 (REACH), annex II.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### General information

In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

#### After inhalation

Remove affected person from the immediate area. Ensure supply of fresh air.

#### After skin contact

Wash off immediately with soap and water.

#### After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.).

#### After ingestion

Never give anything by mouth to an unconscious person. If individual is drowsy or unconscious, place in recovery position (on left side, with head down). Rinse out mouth and give plenty of water to drink.

Trade name: Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution

Current version : 8.0.0, issued: 22.03.2013

Replaced version: 7.1.0, issued: 14.03.2013

Region: GB

**4.2 Most important symptoms and effects, both acute and delayed**

No data available.

**4.3 Indication of any immediate medical attention and special treatment needed**

No data available.

**SECTION 5: Firefighting measures**

**5.1 Extinguishing media**

**Suitable extinguishing media**

Product itself is non-combustible; adapt fire extinguishing measures to surrounding areas.

**Unsuitable extinguishing media**

No data available.

**5.2 Special hazards arising from the substance or mixture**

None known

**5.3 Advice for firefighters**

Use self-contained breathing apparatus. Wear protective clothing.

**SECTION 6: Accidental release measures**

**6.1 Personal precautions, protective equipment and emergency procedures**

**For non-emergency personnel**

Refer to protective measures listed in sections 7 and 8. Ensure adequate ventilation.

**For emergency responders**

No data available. Personal protective equipment (PPE) - see Section 8.

**6.2 Environmental precautions**

Do not discharge into the drains/surface waters/groundwater.

**6.3 Methods and material for containment and cleaning up**

Pick up with absorbent material (e.g., sand, sawdust, general-purpose binder). Send in suitable containers for recovery or disposal.

**6.4 Reference to other sections**

No data available.

**SECTION 7: Handling and storage**

**7.1 Precautions for safe handling**

**Advice on safe handling**

Ensure adequate ventilation.

**General protective and hygiene measures**

Do not eat, drink or smoke during work time. Keep away from foodstuffs and beverages. Avoid contact with eyes and skin. Remove soiled or soaked clothing immediately. Wash hands before breaks and after work.

**Advice on protection against fire and explosion**

No special measures necessary.

**7.2 Conditions for safe storage, including any incompatibilities**

**Technical measures and storage conditions**

Keep container tightly closed in a cool, well-ventilated place.

**Recommended storage temperature**

Value 15 - 25 °C

**Requirements for storage rooms and vessels**

Containers which are opened must be carefully closed and kept upright to prevent leakage. Always keep in containers of same material as the original one.

# Safety data sheet in accordance with 1907/2006/EC

**Trade name:** Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution

**Current version :** 8.0.0, issued: 22.03.2013

**Replaced version:** 7.1.0, issued: 14.03.2013

**Region:** GB

**Advice on storage assembly**

Store the foodstuffs separately.

**7.3 Specific end use(s)**

No data available.

## SECTION 8: Exposure controls/personal protection

**8.1 Control parameters**

No parameters available for monitoring.

**8.2 Exposure controls**

**Appropriate engineering controls**

No data available.

**Personal protective equipment**

**Respiratory protection**

No special measures required.

**Eye / face protection**

Safety glasses (EN 166)

**Hand protection**

Protective gloves (recommended)

**Other**

Standard laboratory clothing.

**Environmental exposure controls**

No data available.

## SECTION 9: Physical and chemical properties

**9.1 Information on basic physical and chemical properties**

<b>Form/Colour</b>			
liquid			
colourless			
<b>Odour</b>			
odourless			
<b>Odour threshold</b>			
No data available			
<b>pH value</b>			
Value	7.24		
Reference temperature	20	°C	
<b>Boiling point / boiling range</b>			
Value	appr.	100	°C
Reference pressure		1013	mbar

# Safety data sheet in accordance with 1907/2006/EC

**Trade name:** Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution

**Current version :** 8.0.0, issued: 22.03.2013

**Replaced version:** 7.1.0, issued: 14.03.2013

**Region:** GB

<b>Melting point / melting range</b>	
No data available	
<b>Decomposition point / decomposition range</b>	
No data available	
<b>Flash point</b>	
No data available	
<b>Auto-ignition temperature</b>	
No data available	
<b>Oxidising properties</b>	
No data available	
<b>Explosive properties</b>	
No data available	
<b>Flammability (solid, gas)</b>	
No data available	
<b>Lower flammability or explosive limits</b>	
No data available	
<b>Upper flammability or explosive limits</b>	
No data available	
<b>Vapour pressure</b>	
No data available	
<b>Vapour density</b>	
No data available	
<b>Evaporation rate</b>	
No data available	
<b>Relative density</b>	
No data available	
<b>Density</b>	
No data available	
<b>Solubility in water</b>	
Remarks	soluble
<b>Solubility(ies)</b>	
No data available	
<b>Partition coefficient: n-octanol/water</b>	
No data available	
<b>Viscosity</b>	
No data available	

## 9.2 Other information

<b>Other information</b>
No data available.

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

Dangerous reactions are not expecting handling the product according to its intended use.

### 10.2 Chemical stability

Stable under recommended storage and handling conditions(See section 7).

# Safety data sheet in accordance with 1907/2006/EC

**Trade name:** Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution

**Current version :** 8.0.0, issued: 22.03.2013

**Replaced version:** 7.1.0, issued: 14.03.2013

**Region:** GB

## 10.3 Possibility of hazardous reactions

None, if handled according to order.

## 10.4 Conditions to avoid

None known

## 10.5 Incompatible materials

none, if handled according to order

## 10.6 Hazardous decomposition products

No hazardous decomposition products if stored and handled as prescribed.

## SECTION 11: Toxicological information

### 11.1 Information on toxicological effects

<b>Acute oral toxicity</b>
No data available
<b>Acute dermal toxicity</b>
No data available
<b>Acute inhalational toxicity</b>
No data available
<b>Skin corrosion/irritation</b>
No data available
<b>Serious eye damage/irritation</b>
No data available
<b>Respiratory or skin sensitisation</b>
No data available
<b>Germ cell mutagenicity</b>
No data available
<b>Reproduction toxicity</b>
No data available
<b>Carcinogenicity</b>
No data available
<b>STOT-single exposure</b>
No data available
<b>STOT-repeated exposure</b>
No data available
<b>Aspiration hazard</b>
No data available
<b>Delayed and immediate effects as well as chronic effects from short and long-term exposure</b>
When handled appropriately adverse health effects are not to be expected.
<b>Other information</b>
The product was classified in toxicological terms on the basis of the results of the calculation procedure outlined within General Directive on Preparations (1999/45/EC).

Trade name: Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution

Current version : 8.0.0, issued: 22.03.2013

Replaced version: 7.1.0, issued: 14.03.2013

Region: GB

## SECTION 12: Ecological information

### 12.1 Toxicity

<b>Fish toxicity</b>
No data available

<b>Daphnia toxicity</b>
No data available

<b>Algae toxicity</b>
No data available

<b>Bacteria toxicity</b>
No data available

### 12.2 Persistence and degradability

No data available.

### 12.3 Bioaccumulative potential

No data available.

### 12.4 Mobility in soil

No data available.

### 12.5 Results of PBT and vPvB assessment

No data available.

### 12.6 Other adverse effects

No data available.

### 12.7 Other information

<b>Other information</b>
Do not discharge into the drains or waters and do not store on public depositories.

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

**Product**

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

**Packaging**

Contaminated packaging should be emptied as far as possible and after appropriate cleansing may be taken for reuse. Packaging that cannot be cleaned should be disposed of in agreement with the regional waste disposal company.

# Safety data sheet in accordance with 1907/2006/EC



**Trade name:** Glucose/Lactate system solution and Glucose/Lactate hemolyzing solution

**Current version :** 8.0.0, issued: 22.03.2013

**Replaced version:** 7.1.0, issued: 14.03.2013

**Region:** GB

## SECTION 14: Transport information

### 14.1 Transport ADR/RID/ADN

The product is not subject to ADR/RID/ADN regulations.

### 14.2 Transport IMDG

The product is not subject to IMDG regulations.

### 14.3 Transport ICAO-TI / IATA

The product is not subject to ICAO-TI / IATA regulations.

### 14.4 Other information

No data available.

### 14.5 Environmental hazards

Information on environmental hazards, if relevant, pls. see 14.1 - 14.3.

### 14.6 Special precautions for user

No data available.

### 14.7 Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code

Not relevant

## SECTION 15: Regulatory information

### 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture EU regulations

#### **Council Directive 96/82/EC on the control of major-accident hazards involving dangerous substances**

Remarks

Annex I, part 1 + 2: not mentioned. With regard to possibly appropriate decomposition products see Chapter 10.

### 15.2 Chemical safety assessment

No data available.

## SECTION 16: Other information

#### **Sources of key data used to compile the data sheet:**

EC Directive 67/548/EC resp. 1999/45/EC as amended in each case.

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

EC Directives 2000/39/EC, 2006/15/EC, 2009/161/EC

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding chapter.

#### **Department issuing safety data sheet**

UMCO Umwelt Consult GmbH

Georg-Wilhelm-Str. 183 , D-21107 Hamburg

Tel.: +49 40 / 79 02 36 300 Fax: +49 40 / 79 02 36 357 e-mail: umco@umco.de

This information is based on our present state of knowledge and experience.

The security data sheet describes products with a view to the security requirements.

However, it should not constitute a guarantee for any specific product properties and shall not establish a legally valid relationship.