



## QuStick™ Strep A

Early diagnosis of Strep A infection means more cost-effective treatment and increased patient satisfaction

- Enhanced sensitivity
- Superior culture classification
- Easy-to-use
- Results within five minutes

# QuStick™ Strep A Rapid Test

Rapid, qualitative detection of group A streptococcal antigen directly from patient throat specimen.

## Enhanced sensitivity

97% sensitivity compared to culture

QuStick™ Strep A's ultra-sensitivity means that it requires fewer colonies to yield a positive result. Which means that even patients presenting less obvious symptoms can be detected

## Superior culture classification (low counts)

91 % @ Rare (< 5 colonies)

100 % @ 1+ (< 10 colonies)

## Easy-to-use

CLIA waived, requires minimal training

## Fast

Results in five minutes

## Economical. Convenient.

QuStick™ eliminates waste with room temperature storage and a two year shelf life from date of manufacture

## Complete

The kit includes everything you need

## Each kit includes:

- Individually pouched test strips
- Reagent A and reagent B
- Positive and negative external controls
- Sterile throat swabs
- Extraction tubes
- Instructions for use
- Handy pictorial quick guide showing step-by-step procedure



Reference Nos.

6000-025 (25 tests)

6000-050 (50 tests)

EKF Diagnostics  
1261 North Main Street  
Boerne, TX 78006

☎ 1 800 531 5535

☎ (830)249 0772

✉ sales2@ekfdiagnostics.com

[www.ekfusa.com](http://www.ekfusa.com)

**CLIA  
WAIVED**

CODE\*:  
87880QW

\* Exemplar code provided. Refer to regional CMS fee schedule for relevant codes and reimbursement rates.