



# EKF Link

Connectivity at the point-of-care

Scan here to view



product webpage

Easy-to-use interface

Safe and secure

Connect multiple point-of-care analyzers

Maximum flexibility

**EKF**  
Diagnostics

# Connect all of your point-of-care analyzers with EKF Link

	Feature	Advantage
<b>Safe and secure</b>	Integrated device connector	Connect multiple Point of Care analyzers securely and easily to LIMS/EMR.
<b>Maximum flexibility</b>	Advanced architecture	Hosted on your laptop, or server. <b>Your data, your choice.</b> Data accessible from most web browser using security protocol without installing any additional module.
<b>Test integrity</b>	Patient and QC tests management	Manage test in table or/and charts Check results Review exceptions Add comments
<b>Traceability</b>	Operator management	Manage Operator certification which reduces the risk for errors and helps to meet strict regulatory and quality requirements.
<b>Visibility</b>	Instrument management	Manage your instrument remotely, check their status, so you know if a device has to be replaced. Remote configuration, so you don't need to go physically on each instrument to change their settings.
<b>Archive</b>	Report management	You can export test results on pdf and csv report for further review.

EKF Link is a powerful middleware solution for healthcare providers, centralizing data management for point-of-care (POC) analyzers on one platform.

A flexible, open data management solution that seamlessly integrates with most vendors' point-of-care analyzers, allowing real-time management of patient test results, quality control (QC) data, operator management, and analyzer configuration.

Designed to provide a flexible and secure method for connecting point-of-care devices to hospital and laboratory IT systems, enabling the transfer of data from the analyzer to the hospital or laboratory database for subsequent processing. As an optional software platform, it can be installed on a PC or a server, allowing authorized users to access it anywhere at any time.

Catalog no.	Name
0760-0102	EKF Link license
0760-0104	EKF Link Lite

## EKF Link Intended Use

EKF Link is a data and device management system that enables healthcare professionals to transfer, manage, store and display measurement results from point-of-care diagnostic testing analyzers. EKF Link enables the allocation of additional information to measurement results and facilitates the export or transfer of data for documentation purposes.

EKF Link is intended to be used as an optional software application with EKF point-of-care analyzers and open to be interfaced to other vendor's point-of-care testing analyzers verified for connection with EKF Link. EKF Link facilitates the setting of configurable options of connected analyzers.

EKF Link is not intended for use in patient diagnostics, screening, monitoring or for making treatment decisions. Such decisions should be made based on the results displayed on the POC device.



## Manufacturer

EKF-diagnostics GmbH  
Ebendorfer Chaussee 3  
39179 Barleben Germany  
Tel: +49 (0) 39203 511 0  
Email: sales@ekfdiagnostics.com  
**ekfdiagnostics.com**

## Distributed by