ekfdiagnostics.com

Powered by EKF Link









Quo-Test®

Easy and reliable HbA1c measurement

Scan here to view product webpage

- Fully automated analyzer
- Results in just four minutes
- Just 4µl blood sample size
- Connection to EKF Link with Quo-Connect

EKF Diagnostics

Quo-Test®

Easy-to-use

- Intuitive, step-by-step on screen instructions minimise training time
- Multilingual menu
- Just 4µl blood required from finger prick or venous sample

Efficient data handling

- Barcode reader included. Only one scan per lot of cartridges required
- Stores up to 7000 results for efficient patient data management
- User-selectable dual reporting. User can select % DCCT, IFCC mmol/mol, eAG mg/dl or eAG mmol/l.
- USB Port
- Optional label printer available
- Powered by EKF Link and Quo-Connect for transfer of results into LIS/LIMS


Fast and accurate

- Results within four minutes
- Boronate affinity method which is recognized as interference free
- Measuring range: 4–15% A1c DCCT
- Imprecision: CV <3% at 7% A1c DCCT
- Unaffected by Hb variants, which do not result in reduced erythrocyte life span
- Traceable to the IFCC reference method

Ergonomic design

- Compact: 205 mm(h) x 205 mm(w) x 135 mm(d)
- Lightweight: 1.3kg

Catalog no.	Name
3121-0011_VS	• Quo-Test® A1c analyzer
6002-0011_VS	• Quo-Connect Start Set
0760-0102	• EKF Link license



Collect 4 µl blood sample.



Insert sample collector into the Quo-Test cartridge.



Insert cartridge into analyser. Results displayed within 4 minutes.

Manufacturer
EKF-diagnostic GmbH
Ebendorfer Chaussee 3
39179 Barleben Germany
Tel: +49 (0) 39203 511 0
Email: sales@ekfdiagnostics.com
ekfdiagnostics.com

Distributed by

QuoTest Data EN EU Rev 4.0 10-2026

This document is for general information purposes only. A product may only be operated or used as per the instructions/manuals provided with it.